

ChemRisk/Shonka Research Associates, Inc., Document Request Form

(This section to be completed by subcontractor requesting document)

J. Lamb / 1034A
Requestor Document Center (is requested to provide the following document)

Date of request 8/18/95 Expected receipt of document 9/15/95

Document number K2-1278 Date of document none

Title and author (if document is unnumbered)

(This section to be completed by Document Center)

Date request received 8/24/95

Date submitted to ADC 8/30/95

Date submitted to HSA Coordinator 8/24/95

(This section to be completed by HSA Coordinator)

Date submitted to CICO 8/30/95

Date received from CICO 9/6/95

Date submitted to ChemRisk/Shonka and DOE 9/6/95

(This section to be completed by ChemRisk/Shonka Research Associates, Inc.)

Date document received _____

Signature _____

B. E. Roche or George

Distillation, Drying and Testing Purity of Hg.

ENT NO.
1278

General:

Most Mercury that comes in the laboratory for treatment, contains substances which effect its use as a reagent therefore it is necessary to follow the outline as here-after described, in order that those foreign substances which effect its actions as a reagent may be removed.

Procedure:

1. Initial Wash

Take 25 lbs of impure or contaminated Hg. Place in 1500 or 2000 cc container and wash by passing Hot H₂O through the Hg for at least 4 hrs. depending upon the amount of organic or other foreign substances present. Note; (In case oils are present, this time should be increased to 6 hrs.). The hot water which is passed through should be at such a rate as to break the surface tension of the Hg. This is best obtained by allowing the water to pass through a glass tube and released under the surface of the Hg near the bottom of the bottle or container.

2. Mercury and Water Separation

After completing #1, decant the water above the Mercury and remove all possible water which is on surface, by the use of a vacuum using a trap to collect the water.

3. Removing amalgams and Acid Neutralization

Take the Mercury from No 2 and pass it through a 30 to 40 cm. tube containing 8% HNO₃. Then through a tube of the same size containing distilled water. These solutions should be changed after passing through 10 lbs of Hg. The rate at which the Hg is passed through these solutions should exceed not 5 lbs per hour.

4. First Drying (Drierite method)

The Hg from No 3 is passed over a tube containing Drierite. The tube should be from 12 to 14 in. in diameter and 30 cm. in length. Place glass wool at top and bottom of Drierite column so as to insure diffusion of Mercury allowing a greater surface to be exposed. The Mercury level should be kept at 3 to 5 cms. above the upper glass wool mat. The Mercury should not pass from Drierite tube faster than 5 lbs. per hour. The Drierite should be changed after passing 30 lbs of Hg through the column. Note: (The Mercury should not contain more than 0.3 of 1% of water by weight. See No. 7 for test.)

5. Distillation of Mercury

The Mercury from No. 4 is placed in the open top reservoir of still. (See figure 106, page 685, Analytical Chemistry. Treadwell and Hall, Vol2, Quantitative, 9th, English Edition.) Make sure that the distance from the Mercury level in the reservoir is 75 cms, lower than the upper level of the Hg in the distillation flask. (Distillation flask should be $\frac{1}{2}$ full.) The center tube, or condensation tube which is to receive the Mercury vapors, should always extend 1 cm. above the Mercury level in the distillation flask. CAUTION: These distances vary according to the barometric pressure and the amount of heat applied to the distillation flask. At least 5 cms. should be allowed from the top of the Hg in the reservoir, allowing for pressure caused by Hg when heat is applied, and also any changes in barometric reading

The lower condensation tube, or center tube, should be at least 82 cms in length, measuring from the bottom of the open top Hg reservoir to the

PLANT RECORDS DEPT.
CENTRAL FILES

REC

FILE 127
X-REF.

PLANT RECORDS

outlet arm of the Hg seal flask, at the lower end of center or condensation tube. (Note: Center or condensation tube should extend to within $\frac{1}{2}$ cm. of the bottom of the Hg seal flask, (Lower flask). To begin distillation, connect vacuum pump to outlet of Hg seal tube (lower receiving flask), and pull the Hg by vacuum to the desired level in the distillation flask, (half full). (1 cm. below the upper end of center or condensation tube). When level is reached add Hg to the open top reservoir bringing the Hg to the calibrated level. Often it is necessary to add small quantities of Hg while vacuum pump is in operation in order to bring Hg to desired level in distillation flask.

Now a small flame is started under distillation flask. (CAUTION: This should be watched constantly, so that the flame is not too high, causing back pressure on the open Hg reservoir, causing the same to overflow.) Under heat the Hg may rise to 1 to 2 cms. in the reservoir but never higher. The vacuum pump remains in operation until the Hg seal flask is filled to within $\frac{1}{2}$ cm. of the outlet arm. Now close the bottom screw clamp, disconnect the vacuum pump, and allow the distilled Hg to rise to not less than 50 cms. in the condensation or center tube. Then open the screw clamp slowly and place receiving vessel beneath outlet arm of Hg sealing flask and allow Hg to continue.

6. Final Drying (H₂SO₄ Method.)

The Hg from 5 is placed on electric hot plate under hold. Connect to vacuum with trap between vacuum and Hg container. Dry air is introduced to Hg container using a glass tube extending to the bottom of the container, in order to prevent bumping. The air is dried through conc. H₂SO₄ with a trap between the acid and the Hg. This process is allowed to continue for 4 hours at 180 to 200 degrees F. Remove and allow to cool and strain through 4 layers of clean gauze.

7. Test for Moisture Content

25 cc. sample of Hg from No. 6 is shaken out with 25 cc. of alcohol which is an aliquote portion of a 500 cc. sample on which a blank has been run, using Fishers reagent. Decant the alcohol, measure and titrate with Fishers reagent. The amount of water present should not exceed 0.01% by weight. (Note: In case Hg contains a higher percentage of moisture return and repeat No. 6.)

8. Qualitative Test for Metal Radicals.

The principal impurities found in Hg are: Copper, cadmium, zinc, and some times silver and gold. These substances are tested for as outlined in Langes Handbook Of Chemistry, 4th addition, Page 946. Should any of the above be present, repeat No. 3 through 7 inclusive. (Note: Increase HNO₃ used in No. 3 to 10%.

9. Labeling and Accounting of Reagent Mercury

The Hg from No. 8 is labeled (Distilled and Dried Hg.) Moisture content, date prepared, and quantity or weight. The Hg is tightly stoppered and stored, and shall be accounted for as prescribed by chief chemist in charge.

J. Geary