

(This section to be completed by subcontractor requesting document)Requestor J. Lamb / 1034A Document Center (is requested to provide the following document)Date of request ~~2/2/96~~ 3/29/96 Expected receipt of document ~~3/16/96~~ 4/30/96Document number K2-1311 Date of document 8/27/45Title and author (if document is unnumbered)

(This section to be completed by Document Center)Date request received 4/8/96Date submitted to ADC 4/15/96Date submitted to HSA Coordinator 4/8/96**(This section to be completed by HSA Coordinator)**Date submitted to CICO 4/15/96Date received from CICO 4/19/96Date submitted to ChemRisk/Shonka and DOE 4/23/96**(This section to be completed by ChemRisk/Shonka Research Associates, Inc.)**

Date document received _____

Signature _____

2748

8-2-45

LABORATORY B

REPORT NO.

KZ 1311

Report on
Tower Discharge Analysis

August 27, 1945

By J. A. Klacsmann

Copies to:

Dr. D. E. Hull ✓
Dr. F. W. Hurd
Mr. C. B. Slade
Mr. K. W. Bahler
Dr. C. M. Birdsall
Mr. J. A. Klacsmann

RECORD COPY
PLANT RECORDS 1950

This document has been approved for release
to the public by: *W. Kelly*

for Arthur S. Zwick
Technical Information Officer
Oak Ridge K-25 Site

Date
4/19/96

Carbide and Carbon Chemicals
Corporation, Operating Contractor for
the U.S. Atomic Energy Commission.

Tower discharge samples were to be analyzed to estimate whether salvage would be economically advisable. Only the order of magnitude of material was needed.

The material came in a 3% basic Sodium carbonate solution, containing approximately 10 mgs of material per liter.

The various methods tried were peroxide precipitation directly, sulfide precipitation directly, and an ammonia precipitation method (See Procedure I). These were each followed by a residue count. Besides these methods plating from the carbonate solution was tried, and the resultant film was counted.

The sulfide and peroxide methods gave unsatisfactory results, and further mention will not be given (See Table I).

The ammonia method worked satisfactorily, and a procedure was written (Procedure I).

Control samples were run as follows: A known amount of material of known counting rate was added to an unknown solution and the extra counts were determined experimentally. It was shown that 40% of the counts added were obtained.

Sample 1542 was run with 550 counts added to it as a control. 10 ml portion of 500 ml volume of 463 mg. of I_2^{131} (as nitrate) was diluted to 100 ml and 10 ml aliquot of that was added to sample (100 ml).

As can be seen from the data only 220 of the 550 counts added were obtained.

The answers were therefore multiplied by $2\frac{1}{2}$ before reporting (3 obtained by use of 1645).

The plating procedure (Procedure II) was worded out, and gave results comparable to the calculated Ammonia method.

Summary:

1. The results will be good the the desired certainty. (See Table I).
2. The plating method is a fast simple method which can be run by routine people.

TOWER DISCHARGE SAMPLES

Report on methods:

TABLE I

Sample No.	Method	Sample W and Mgs.	Volume ml.	Counts Per 100 ml.
1542	N H ₃	3.7	100	133
1542 + 550 $\frac{\text{Counts}}{\text{min}}$	N H ₃	5.5	100	350
1564	N H ₃	2.3	100	95
1564 (Clear)	H ₂ S	49.9	100	11
1564	Plating	—	30	240
1585	N H ₃	2.3	100	62
1585	N ₂ S	49.6	100	17
1594	N H ₃	5.6	100	64
1594	Plating	—	100	235
1645	N H ₃	5.9	100	64
1645 + 550 $\frac{\text{Counts}}{\text{min}}$	N H ₃	4.4	100	330
1645	Plating(HN ^o ₃)	—	100	35
1786	Plating	—	100	73
1786	Plating	—	100	71
1786	Plating	—	100	64
1748	Plating	—	100	79
1748	Plating	—	100	67
1748	Plating	—	100	81

69

76

PROCEDURE I

PROCEDURE FOR NH_4OH PRECIPITATION ON DISPOSAL PLANT TOWER DISCHARGE

- I. Acidify 100 ml of solution with nitric acid until solution shows acid.
- II. Evaporate acidified solution to dryness.
- III. Add 10 ml. concentrated nitric acid, washing down sides of beaker with acid.
- IV. Evaporate until dry and bake until residue stops omitting acid fumes.
- V. Dissolve the residue in a minimum of CO_2 free water.
- VI. Heat solution just to boiling and add an equal volume of CO_2 free 4 M NH_4OH .
- VII. Centrifuge for 30 seconds.
- VIII. Pour off excess liquid.
- IX. Add concentrated nitric acid dropwise until precipitate dissolves being careful not to use an excess of acid.
- X. Add 4 M NH_4OH in amount equal to solution of sample.
- XI. Centrifuge for 30 seconds and pour off excess liquid.
- XII. Dry precipitate in centrifuge tube under infrared - radiator until the ppt. becomes hard.
- XIII. Loosen the dried precipitate with a platinum spotula and transfer to a weighted platinum dish.
- XIV. Ignite sample for 15 minutes and weigh.
- XV. Dissolve ignited sample in nitric acid.
- XVI. Evaporate solution to dryness.
- XVII. Heat sample gently until it turns to a red color.
- XVIII. Add 4 M HCL to dissolve.
- XIX. Transfer sample to silver foil disc for correcting as follows:
 - a. Place disc in bottom of petri dish placed on a hot plate.
 - b. Drop sample solution on disc from a dropper no faster than the solution will evaporate to form a deposit in the middle portion of the disc.
 - c. Rinse platinum dish twice with some 4 M HCL.
 - d. Ignite disc in oven for one minute at 25°C .